

EXPERIMENT 5  
ACID AND BASE STRENGTH

**PURPOSE:**

1. To distinguish between acids, bases and neutral substances, by observing their effect on some common indicators.
1. To distinguish between strong and weak acids and bases, by conductivity testing.
3. To identify an unknown, as an acid (strong or weak), a base (strong or weak) or a neutral substance.

**PRINCIPLES:**

We frequently encounter acids and bases in our daily life.

**Acids** were first associated with the sour taste of citrus fruits. In fact, the word **acid** comes from the Latin word **acidus**, which means, “**sour**”. Vinegar tastes sour because it is a dilute solution (about 5 percent) of acetic acid.

Citric acid is responsible for the sour taste of a lemon. Hydrochloric acid is the acid in the in the gastric fluid in your stomach, where it is secreted at a strength of about 5 percent.

A normal diet provides mostly acid-producing foods. Water solutions of acids are called **Acidic Solutions**.

**Bases** have usually a **bitter** taste and a slippery feel, like wet soap. The bitter taste of tonic water comes from a natural base, quinine. Common medicinal antacides (used to relieve heartburn) and bitter tasting Milk of Magnesia, a common laxative, (a suspension of about 8 percent of magnesium hydroxide) are bases. The most important of the strong bases is sodium hydroxide, a solid whose aqueous solutions are used in the manufacture of glass and soap.

Water solutions of bases are called alkaline solutions or **Basic Solutions**.

Substances used to determine whether a solution is acidic or basic are known as **indicators**.

Indicators are organic compounds that change color in a specific way, depending on the acidic or basic nature of the solution. A wide variety of indicators are commonly used in the chemistry laboratory, to identify the acidic or basic nature of an aqueous solution. This experiment uses only two types of indicators: **Litmus** (a vegetable dye) and **Phenolphthalein**.


Some of the characteristic properties commonly associated with acids and bases in aqueous solutions are the following:

ACIDS	BASES
Sour taste	Bitter taste
Change the color of Blue Litmus Paper To <b>RED</b>	Change the color of Red Litmus Paper to <b>BLUE</b>
Do not change the color of phenolphthalein	Change the color of phenolphthalein to <b>PINK</b>
React with carbonates to produce CO <sub>2</sub>	React with acids
React with bases	React with acids

When acids and bases react with one another in equal proportions, the result is a **neutralization reaction**, which produces neutral products: salt and water. Neutral means in this context, that these products do not change the color of litmus or phenolphthalein, do not have a sour or bitter taste, therefore they are “neither acidic nor basic”.

EXPERIMENT 5  
ACID AND BASE STRENGTH

The following equation represents a typical acid-base neutralization reaction:



**A salt is any compound of a cation (other than  $\text{H}^+$ ) with an anion (other than  $\text{OH}^-$  or  $\text{O}^{2-}$ ).**

Acid properties are often opposite to base properties, and vice versa; a base is an anti-acid, and an acid is an anti-base.


Several theories have been proposed to answer the question “What is an Acid or a Base”?

One of the earliest and most significant of these theories was proposed by a Swedish scientist, Svante Arrhenius in 1884.

**According to Arrhenius:**

AN ACID	A BASE
<b>Is a hydrogen-containing substance that dissociates to produce hydrogen ions, <math>\text{H}^+</math>, in aqueous solutions.</b>	<b>Is a hydroxide-containing substance that dissociates to produce hydroxide ions, <math>\text{OH}^-</math>, in aqueous solutions.</b>
<b>The hydrogen ions, <math>\text{H}^+</math>, are produced by the dissociation of acids in water</b> $\text{HA(aq)} \longrightarrow \text{H}^+(\text{aq}) + \text{A}^-(\text{aq})$ <b>Acid</b>	<b>The hydroxide ions, <math>\text{OH}^-</math>, are produced by the dissociation of bases in water</b> $\text{MOH(aq)} \longrightarrow \text{M}^+(\text{aq}) + \text{OH}^-(\text{aq})$ <b>Base</b>
<b>An ACID SOLUTION contains an excess hydrogen ions, <math>\text{H}^+</math>.</b>	<b>A BASE SOLUTION contain an excess of hydroxide ions, <math>\text{OH}^-</math>.</b>
<b>Examples: <math>\text{HCl(aq)}</math>, <math>\text{HC}_2\text{H}_3\text{O}_2(\text{aq})</math></b>	<b>Examples: <math>\text{NaOH(aq)}</math>, <math>\text{NH}_4\text{OH(aq)}</math></b>

Today we know that  $\text{H}^+$  ions cannot exist in water, because a  $\text{H}^+$  ion is a bare proton, and a charge of +1 is too concentrated for such a tiny particle. Because of this, any  $\text{H}^+$  ion in water immediately combines with a  $\text{H}_2\text{O}$  molecule to form a hydrated hydrogen ion,  $\text{H}_3\text{O}^+$  [that is,  $\text{H}(\text{H}_2\text{O})^+$ ], commonly called a **hydronium** ion.



While it is a known fact that the hydrogen ion does not exist alone, as  $\text{H}^+(\text{aq})$ , but is stable in aqueous solution in the form of the hydronium ion,  $\text{H}_3\text{O}^+(\text{aq})$ , it is an accepted simplification to represent the hydronium ion,  $\text{H}_3\text{O}^+(\text{aq})$ , as a hydrogen ion,  $\text{H}^+(\text{aq})$


EXPERIMENT 5  
ACID AND BASE STRENGTH

In beginning courses, formulas for acids (and no other compounds except water) are written with the dissociable hydrogen atoms (**acidic hydrogen atoms**) first, as in:

**H**Cl (aq)      Hydrochloric acid

**H**C<sub>2</sub>H<sub>3</sub>O<sub>2</sub>(aq) Acetic acid

Only the hydrogen written first is capable of being released as hydrogen ion, H<sup>+</sup>; the other three hydrogen atoms do not yield H<sup>+</sup> ions in aqueous solution

Methane	CH <sub>4</sub>		Are not acids, since do not provide H <sup>+</sup> ions to aqueous solutions. Their hydrogen atoms are therefore not written first in their formulas.
Ammonia	NH <sub>3</sub>		
Urea	NH <sub>2</sub> —CO—NH <sub>2</sub>		
Glucose	C <sub>6</sub> H <sub>12</sub> O <sub>6</sub>		

With a slight modification (the introduction of the H<sub>3</sub>O<sup>+</sup> ion), the Arrhenius definitions of acid and base are still valid today, as long as we are referring to aqueous solutions.


**Arrhenius definitions of ACIDS and BASES adapted to aqueous solutions**

When we dissolve <b>an acid (a molecular substance)</b> in water, the molecules of acid react with water to produce <b>H<sub>3</sub>O<sup>+</sup> ions</b> .	When we dissolve <b>a base (an ionic substance)</b> in water, the metallic ion and the hydroxide ions, OH <sup>-</sup> separate.
H <sub>2</sub> O	H <sub>2</sub> O
HCl(g) $\xrightarrow{\text{H}_2\text{O}}$ H <sub>3</sub> O <sup>+</sup> (aq) + Cl <sup>-</sup> (aq) Accepted simplification: HCl(aq) $\xrightarrow{\text{H}_2\text{O}}$ H <sup>+</sup> (aq) + Cl <sup>-</sup> (aq)	NaOH(s) $\xrightarrow{\text{H}_2\text{O}}$ Na <sup>+</sup> (aq) + OH <sup>-</sup> (aq) Accepted simplification: NaOH(aq) $\xrightarrow{\text{H}_2\text{O}}$ Na <sup>+</sup> (aq) + OH <sup>-</sup> (aq)
The aqueous solution of hydrochloric acid contains ions only (no molecules)	The aqueous solution of sodium hydroxide contains ions only (no molecules)
The acidic solution is a strong electrolyte. (Complete dissociation took place).	The basic solution is a strong electrolyte. (Complete dissociation took place)
<b>Acids, completely dissociated in ions</b> in aqueous solutions, are called <b>strong acids</b> .	Soluble metallic hydroxides, <b>completely separated</b> in aqueous solution, are called <b>strong bases</b> .

Other substances, although they do in fact produce hydrogen ions, H<sup>+</sup>, when dissolved in water, dissociate only partially.

**Substances that produce H<sup>+</sup> ions by partial dissociation are called weak acids. It follows that weak acids are weak electrolytes.**

For example, acetic acid, HC<sub>2</sub>H<sub>3</sub>O<sub>2</sub>, found in vinegar, is a weak acid.



The double arrows in the above equation for the partial dissociation of acetic acid indicate that the dissociation reaction for this substance reaches **equilibrium**.



**EXPERIMENT 5**  
**ACID AND BASE STRENGTH**

In 1923 new definitions of acids and bases were proposed simultaneously by Bronsted and Lowry. The Bronsted/Lowry theory of acids and bases extends the Arrhenius definitions to more general situations, which explain the behavior of weak bases and do not require the solvent to be water.

**According to the Bronsted/Lowry theory:**


**AN ACID:**

**Is a proton,  $H^+$ , donor**

**A BASE:**

**Is a proton,  $H^+$ , acceptor**

**AN ACID-BASE REACTION (NEUTRALIZATION REACTION) IS THE TRANSFER OF A  $H^+$  ION.**



In summary, both acids and bases have characteristic properties and can either be strong or weak, as shown below:

	<b>ACIDS</b>		<b>BASES</b>	
Arrhenius definition	produce $H^+$ ions in aqueous solution		produce $OH^-$ ions in aqueous solution	
Bronstead/Lowry definition	$H^+$ donors		$H^+$ acceptors	
Electrolyte Strength	STRONG ACIDS (strong electrolytes)	WEAK ACIDS (weak electrolytes)	STRONG BASES (strong electrolytes)	WEAK BASES (weak electrolytes)
Extent of dissociation	completely dissociated	partially dissociated	completely dissociated	partially dissociated
Symbols used to show extent of dissociation	$\longrightarrow$	$\rightleftharpoons$	$\longrightarrow$	$\rightleftharpoons$
Particles present solution	ions only	mostly molecules and a few ions	ions only	mostly molecules and a few ions

---

EXPERIMENT 5  
ACID AND BASE STRENGTH

Keep in mind that **strong** and **concentrated** are not interchangeable terms when applied to acids and bases:

**STRONG:** refers to the extent to which an acid or base dissociates in water.

**CONCENTRATION:** describes how much of an acidic or basic compound is present in a solution.

Strong acids and bases are 100% dissociated. Therefore we cannot interpret the relative acidic or basic strength among strong acids and bases.

**Strong Acids and Strong Bases have acidic and basic properties to an extreme.**

The situation is quite different for weak acids and bases. They dissociate partially and to different degrees. The more an acid dissociates, the more free  $H^+$  ions are present and the stronger the weak acid.

A similar situation exists for weak bases: the more a weak base dissociates, the more free  $OH^-$  ions are present and the stronger the weak base. No attempt is made in this experiment to rank the weak acids or bases according to their relative strength.

In this experiment:

- Indicators will be used to distinguish between acids, bases and neutral substances, and
- A semi-quantitative conductance testing will be performed to distinguish between strong and weak acids. The result of the conductance test should clearly distinguish between strong and weak electrolytes, and will identify the acid or the base as strong or weak.

### **PROCEDURE:**

In this experiment, the effect on indicators and the conductance of six aqueous solutions and an unknown will be determined. The unknown, assigned individually, is one of the six solutions tested. Two types of indicators will be used to determine the acidic or basic character of the solutions: **Litmus paper** (Red and Blue) and **Phenolphthalein**.

From the data you gather, you will be able to determine: the acidic, basic, or neutral character of the solution, the electrolyte character and the formula of the predominant species in solution. If the solution is acidic or basic, you will be able to determine if the acid or the base is strong or weak. All aqueous solutions tested (including the unknown) have the same concentration: 0.10 M. The formulas and the names of the six known solutions tested are listed below:


1.  $HC_2H_3O_2(aq)$       **acetic acid**
2.  $NH_3(aq)$             **aqueous ammonia**
3. **D.I.  $H_2O$**             **deionized water**
4.  $HCl(aq)$              **hydrochloric acid**
5.  $NaOH(aq)$           **sodium hydroxide**
6.  $NaCl(aq)$             **sodium chloride**

---

EXPERIMENT 5  
ACID AND BASE STRENGTH

**I. Indicator testing****1. Litmus paper Test**

- (a) Lay a sheet of regular white paper or a rectangular clean and dry sheet of paper towel on your working area.
- (b) Obtain 8 pieces of Red Litmus paper and 8 pieces of Blue Litmus paper.  
When using litmus paper (Red or Blue), there is no need to use a full strip of indicator paper.  
Avoid wasting indicator paper by using only one-half or one-third of indicator paper strip for each test.  
Note that:
- The color of “RED” LITMUS PAPER is actually “PINK”
  - The color of “BLUE” LITMUS PAPER is actually “FAINT LAVENDER”.
- (c) Arrange and label the pieces of litmus paper on the sheet of white paper.



EXPERIMENT 5  
ACID AND BASE STRENGTH



EXPERIMENT 5  
ACID AND BASE STRENGTH

- (d) Apply 1-2 drops of each solution (6 known solutions and two unknown solutions) on each strip of litmus paper.
- If the solutions are provided in dropper bottles:
    - Dispense 5 solutions directly from the dropper bottles, and
    - Dispense the D.I. H<sub>2</sub>O and the two unknown solutions by using disposable plastic droppers.
  - If the solutions are NOT provided in dropper bottles:
    - You may dispense the required 1-2 drops of solutions onto the litmus paper strips by:
      - Using disposable plastic droppers. You may not use the same disposable plastic dropper when switching for one solution to another.
      - OR
      - Using the narrow hollow glass tube from your locker.  
The narrow hollow glass tube must be washed thoroughly when switching from one solution to another.
- (e) Record your observations in the Laboratory Notebook.

## 2. Phenolphthalein Test

(a) Wash (tap water followed by D.I water) nine test tubes and allow the water to drain.

(b) Label the test tubes as follows:

<b>1</b>	<b>2</b>	<b>3</b>	<b>4</b>	<b>5</b>	<b>6</b>	<b>7</b>	<b>8</b>
HC <sub>2</sub> H <sub>3</sub> O <sub>2</sub> (aq)	NH <sub>3</sub> (aq)	D.I. H <sub>2</sub> O	HCl(aq)	NaOH(aq)	NaCl(aq)	Unk. # _____	Unk. # _____

- (c) Fill **test tube # 3** with approximately 1 mL of D.I. H<sub>2</sub>O.  
The level of D.I. H<sub>2</sub>O in this test tube may be used to match the level of the remaining seven solutions in each test tube.
- You may use your 10 mL graduated cylinder to measure this volume of D.I. water, OR
  - You may use a disposable plastic dropper to dispense about 20 drops of D.I. water (this volume is approximately equivalent to 1 mL of volume).
- (d) Fill all the other test tubes with the respective solutions, by matching the level of the contents with the level of D.I. H<sub>2</sub>O in test tube # 3.
- If the solutions are provided in dropper bottles dispense the remaining 5 solutions directly from the respective dropper bottles, and dispense the two unknown solutions by using disposable plastic droppers.
  - If the solutions are NOT provided in dropper bottles you may deliver the required volume of each solution (1 mL or about 20 drops) in each test tube by using disposable plastic droppers or a hollow glass tube.
- (e) Add 2 drops of phenolphthalein solution to each of the eight solutions to be tested (Six known solutions and two unknown solutions).
- (f) Observe the color of the solutions and record your observations.


---

EXPERIMENT 5  
ACID AND BASE STRENGTH

## II. Conductivity Testing

### 1. Connecting your “Conductivity Tester Combo”

- (a) Check out a “Conductivity Tester Combo”, composed of a “Power Converter” connected to a “Conductivity Indicator”.
- (b) Make sure that the Power Converter is in the “OFF” position.
- (c) If the dial is not already set to 9V, set the dial to 9 V.
- (d) Keep the dial in this position throughout your experiment.
- (e) Connect the cords:
  - Snap the alligator end of the red cord onto the (+) post (labeled “R”) of the conductivity indicator and the black cord onto the (-) post (labeled “B”) of the conductivity indicator.
  - If not already taped, solidly tape the connections of the cords to the conductivity indicator with electrical tape (available at the instructor’s desk)
  - Insert the other two ends of the cords into the outlets of the Power Converter:
    - Red end of the cord into the red outlet of the Power Converter
    - Black end of the cord into the outlet of the Power Converter
- (f) After ensuring one more time that the Power Converter is in the “OFF” position, plug in the Power Converter into the power supply available at your station.



### 2. Checking the Conductivity Indicator

- (a) Turn the Power Converter “ON”
- (b) Test the circuit by touching the tips of the two copper electrodes to a coin or a piece of metal.
  - If the conductivity indicator is working properly, the LED will light up brightly or blink.
  - If the LED light does not light up, notify your laboratory instructor.
- (c) Check if the copper electrodes are clean by immersing them into D.I. water. (your 250 mL beaker filled with about 150 mL D.I. water may be used for this purpose)
  - The electrodes are clean if the LED does not light up when the electrodes are immersed in D.I water.
  - If the LED lights up, indicating that the electrodes are dirty, discard the D.I water, and immerse the electrodes again in fresh D.I. water.
  - If necessary, rinse the electrodes repeatedly with fresh portions of D.I water, until you get no response from the LED.

**EXPERIMENT 5**  
**ACID AND BASE STRENGTH**

**3. Preparing your Chemplate**

- (a) Wash (tap water followed by D.I water) the Chemplate and allow the water to drain.  
 (b) Plan which solutions you will place in each depression and record this in your Lab Notebook. Note that the depressions are numbered.



**4. Testing the conductivity of the solutions**

- (a) Fill completely the depressions with the eight solutions: Six known solutions and two unknown solutions (one unknown/student)

It takes about 30 drops to completely fill a depression, without overflowing.

- (b) Test the first solution, by immersing the copper electrodes in the solution.

- In order to obtain consistent observations:
  - The position of the conductivity indicator should be perfectly vertical and perpendicular to the Chemplate, and
  - The two electrodes should always be immersed to the same depth.



- The response of the LED will indicate the following types of conductivity:
  - Bright Light or Blinking Light indicates Strong Conductivity (+)
  - Faint Light indicates Weak Conductivity (+/-)
  - No Light indicates No Conductivity (-)
 Observe and record the response of the LED

- (c) Proceed and test the other solutions

- Before testing the next aqueous solution, rinse the two copper electrodes by immersing them in the beaker containing the D.I water.  
 Note that the electrodes are clean if the LED does not light up when the electrodes are immersed in D.I. water.
- Blot the excess water from the electrodes with tissue paper.
- Continue to the next depression and repeat the previous steps of the procedure.

- (d) Interpret your observations according to the guidelines given in the table below:

Conductance	Electrolyte Character	Predominant Species
+	<b>Strong Electrolyte (SE)</b>	<b>Ions</b>
(+/-)	<b>Weak Electrolyte (WE)</b>	<b>Molecules</b>
(-)	<b>Non-Electrolyte (NE)</b>	<b>Molecules</b>

EXPERIMENT 5  
ACID AND BASE STRENGTH

**III. Cleaning up:**

1. Discard your test solutions in the appropriately labeled waste container.
2. Wash your Chemplate with plenty of tap water.
3. Rinse your Chemplate with D.I. water from your wash bottle.
4. Blot the excess water from the Chemplate with clean paper towel
5. Store the Chemplate in your locker.
6. Discard the used strips of litmus paper in the trash can.

EXPERIMENT 5  
ACID AND BASE STRENGTH  
REPORT FORM

NAME: \_\_\_\_\_ Date: \_\_\_\_\_ Partner: \_\_\_\_\_

	Color with Red Litmus Paper*	Color with Blue Litmus Paper**	Color with Phenolphthalein solution	Conductance (+, +/-, or -)	Electrolyte Character (SE, WE, or NE)	Formula of predominant particles	Acid, Base or Neutral***	Strong or Weak
HC <sub>2</sub> H <sub>3</sub> O <sub>2</sub> (aq) 0.10 M								
NH <sub>3</sub> (aq) 0.10 M								
D.I. H <sub>2</sub> O								
HCl 0.10 M								
NaOH(aq) 0.10 M								
NaCl(aq) 0.10 M								
Unknown # <span style="background-color: yellow; display: inline-block; width: 100px; height: 15px;"></span>								

\* The original color of "Red" Litmus paper is actually Pink

- If its color does not change, report: N.C. (No Change)
- If its color changes to Faint Lavender, report "Blue"

\*\* The original color of "Blue" Litmus paper is actually Faint Lavender

- If its color does not change, report: N.C. (No Change)
- If its color changes to Pink, report "Red"

\*\*\* If the solution is neutral, do not complete the last column (Strong or Weak)

EXPERIMENT 5  
ACID AND BASE STRENGTH

EXPERIMENT 5  
ACID AND BASE STRENGTH  
**PART I: A STUDY OF ACIDIC BEHAVIOR**

1. List below the **formulas** and the **names** of the two acids used in this experiment  
Do not forget to include the state designation "aq", after the formula.

FORMULAS                      NAMES

2. Which acid, listed in number (1) above is a **strong acid**?  
Give its formula below:

\_\_\_\_\_

For the strong acid listed above, write an equation that illustrates its  
Ionization/Dissociation reaction:

3. Which acid, listed in number (1) above is a **weak acid**?  
Give its formula below:

\_\_\_\_\_

For the weak acid listed above, write an equation, that illustrates its  
Ionization/Dissociation reaction:

4. What is the essential difference between STRONG ACIDS and WEAK ACIDS, in  
terms of the **particles** they contain in aqueous solution ?

5. What causes acids to behave the same way toward the indicators used in this  
experiment?

EXPERIMENT 5  
ACID AND BASE STRENGTH

**PART II: A STUDY OF BASIC BEHAVIOR**

1. List below the **formulas** and the **names** of the two bases used in this experiment  
Do not forget to include the state designation "aq", after the formula.

FORMULAS                      NAMES

2. Which base, listed in number (1) above is a **strong base**?  
Give its formula below:

\_\_\_\_\_

For the strong base listed above, write an equation that illustrates its  
ionization/dissociation reaction:

3. Which base, listed in number (1) above is a **weak base**?  
Give its formula below:

\_\_\_\_\_

For the weak base listed above, write an equation, that illustrates its  
Ionization/dissociation reaction:

4. What is the essential difference between STRONG BASES and WEAK BASES, in  
terms of the **particles** they contain in aqueous solution?

5. What causes bases to behave the same way toward the indicators used in this  
experiment?



EXPERIMENT 5

ACID AND BASE STRENGTH

**PART III: ASSESSING THE TESTS PERFORMED**

1. You are testing the acid/basic/neutral character of a solution by using blue litmus paper. Assume that only blue litmus paper is available in the laboratory. Upon adding 1 drop of the solution to a strip of blue litmus paper you notice that the color of the blue litmus paper did not change. What do you conclude about the acidic/basic/neutral character of this solution?
  
2. You are testing the acid/basic/neutral character of **another solution** by using red litmus paper. Assume that only red litmus paper is available in the laboratory. Upon adding 1 drop of the solution to a strip of red litmus paper you notice that the color of the red litmus paper did not change. What do you conclude about the acidic/basic/neutral character of this solution?

**PART IV: IDENTIFICATION OF THE UNKNOWN**

1. Is your unknown # \_\_\_\_\_ an **ACID**, a **BASE**, or a **NEUTRAL SUBSTANCE** ?  
\_\_\_\_\_
  
2. (A) If your unknown is an **ACID** or a **BASE**, is it **STRONG** or **WEAK**?  
\_\_\_\_\_
  
- (B) If your unknown is a **NEUTRAL** substance, identify the substance:  
\_\_\_\_\_ (formula)

EXPERIMENT 5  
ACID AND BASE STRENGTH